As per Infection Control, and in coordination with the OR Radiology Supervisor, the protocol for cleaning lead worn in the OR during AGMPs is as follows:

\*\*Please clean and disinfect the lead aprons with AHP wipes (ie. Accel) and hang to dry\*\*

## Doffing & Cleaning OR Lead:

- 1. Remove gloves
- 2. Hand hygiene
- 3. Remove gown
- 4. Hand hygiene
- 5. Don clean gloves
- 6. First wipe down:
  - i) Wipe down lead front
  - ii) 2<sup>nd</sup> person needed to wipe down lead back
- 7. Remove gloves
- 8. Hand hygiene
- 9. Remove lead apron and place on clean table/trolley
- 10. Hand hygiene
- 11. Keep face/head PPE on
- 12. Don new yellow gown & gloves
- 13. Complete 2<sup>nd</sup> wipe down
- 14. Hang on hangers